ULTIMATE GUIDE TO HEARING AIDS PRICES



Types of Hearing Aids



Invisible Hearing Aid







Behind the Ear: Battery Operated



Completely in the Canal

Receiver in Canal:

Battery Operated

Behind the Ear:

Rechargeable



In The Canal: Battery Operated



Receiver in Canal: Rechargeable



Rechargeable Earbuds

PHNAK









Frequently Asked Questions:

- 1. What are the major brands of hearing aids? Signia | Phonak | Widex | ReSound | Oticon | Starkey
- What is the standard warranty on hearing aids?
 2-years warranty. Some premium hearing aids carry 3-4 years of warranty.
- What is the average life span of a hearing aid? Users typically choose to upgrade or move to newer technology between 4-5 years. However, hearing aids can last longer depending on how well you look after them.
- 4. Where do I get started?

Your first step should be to reach out to a hearing care centre with a qualified Audiologist and book an appointment for a hearing assessment (PTA).

ZERO cost EMI options A Home visit facility



9811227269 Scare@centreforhearing.org





1. INVISIBLE HEARING AID (IIC)

- An IIC hearing aid is the most discreet type of Hearing Aid.
- It fits deeper into the ear canal though it only fits certain canal sizes.
- It is best suited to those with moderate hearing loss.








3. IN THE CANAL (ITC)

A. BATTERY OPERATED

- An ITC hearing aid is slightly larger than a CIC hearing aid.
- It has low feedback, like a CIC hearing aid, but can accommodate more features.
- It is suited for mild to severe hearing loss.
- Requires batteries to be replaced from time to time.





B. RECHARGEABLE

- Don't look like hearing aids, more like Apple AirPods.
- No battery changes, no additional running costs.
- 3 hours full charge, gives 15-18 hours of battery life.
- Requires charging, similar to a mobile phone.





4. RECEIVER IN CANAL(RIC)

A. BATTERY OPERATED


- An RIC Hearing Aid, is similar to a BTE (Behind The Ear) hearing aid, but with the receiver placed within the ear canal.
- It is lighter and prevents the 'plugged-up feeling' or the occlusion effect.
- It is almost invisible.
- It is suitable for mild to severe hearing losses.
- Requires batteries to be replaced from time to time.





B. RECHARGEABLE

- No battery changes, no additional running costs.
- 3 hours full charge, gives 15-18 hours of battery life.
- Requires charging, similar to a mobile phone.






5. BEHIND THE EAR (BTE)

A. BATTERY OPERATED


- Useful for all types of hearing loss.
- It is capable of providing more amplification as compared to other types.
- It is the best choice for hearing impaired children.
- It is larger in size hence more visible, but easier to handle, is sturdy and has a longer battery life.
- Requires batteries to be replaced from time to time.





B. RECHARGEABLE

- Suitable for all types of hearing loss.
- Long battery life for continuous use
- Fast charging.
- User-friendly design with easy controls and smartphone connectivity.
- Enhanced performance with noise reduction and wireless streaming.
- No battery replacements are required.





6. RECHARGEABLE EARBUDS

- Stylish Hearing Aids that don't look like Hearing Aids.
- Similar looking to Apple Air-pods or similar Bluetooth devices.
- Suitable for mild to severe hearing loss.
- Easy rechargeability with upto 26 hours use per charge.
- Portable charging case.
- Direct streaming from iPhone and iPad or android smartphones.



For Competitive Prices and **Special Offers** on **Hearing Aids**

Contact Centre for Hearing

www.centreforhearing.org

9811227269

Financing Options

Centre for Hearing has associations with financing companies offering 6-9-11-18 month 0% EMIs.