Ultimate Guide to Hearing Aid Prices in India | 2025



Types of Hearing Aids



Invisible Hearing Aid



Completely in the Canal



In The Canal: **Battery Operated**



In Thehe Canal: Recharargeable



Receiver in Canal: **Battery Operated**



Receiver in Canal: Rechargeable





Behind the Ear: Rechargeable



Rechargeable **Earbuds**














Freque ently Asked Questions:

ire the major brands of hearing aids?

- What a Phonak | Widex | ReSound | Oticon | Starkey Signia |
 - s the standard warranty on hearing aids?
- 2. What is warranty. Some premium hearing aids carry 3-4 years of warranty. 2-years
 - s the average life span of a hearing aid?
- What is ypically choose to upgrade or move to newer technology between Users tyrs. However, hearing aids can last longer depending on how well you 4-5 yecter them. look aft
 - do I get started?
- Where 'st step should be to reach out to a hearing care centre with a Your fired Audiologist and book an appointment for a hearing assessment qualifie









1. INVISIBLE HEARING AID (IIC)

- An IIC hearing aid is the most discreet type of Hearing Aid.
- It fits deeper into the ear canal though it only fits certain canal sizes.
- It is best suited to those with moderate hearing loss.
- Now available as a ready-to-fit, rechargeable hearing aid





2. COMPLETELY IN THE CANAL (CIC)

- A CIC hearing aid is also discreet and fits into the ear canal.
- It is easy to use for telephonic conversations.
- It is quite effective in blocking out background noise.











3. IN THE CANAL (ITC)

A. BATTERY OPERATED

- An ITC hearing aid is slightly larger than a CIC hearing aid.
- It has low feedback, like a CIC hearing aid, but can accommodate more features.
- It is suited for mild to severe hearing loss.
- Requires batteries to be replaced from time to time.





B. RECHARGEABLE

- Don't look like hearing aids, more like Apple AirPods.
- No battery changes, no additional running costs.
- 3 hours full charge, gives 15-18 hours of battery life.
- Requires charging, similar to a mobile phone.















4. RECEIVER IN CANAL (RIC)

A. BATTERY OPERATED



- It is lighter and prevents the 'plugged-up' feeling' or the occlusion effect.
- It is almost invisible.
- It is suitable for mild to severe hearing losses.
- Requires batteries to be replaced from time to time.






Popular

B. RECHARGEABLE

- No battery changes, no additional running costs.
- 3 hours full charge, gives 15-18 hours of battery life.
- Requires charging, similar to a mobile phone.














5. Behind the Ear (BTE)



A. BATTERY OPERATED

- Useful for all types of hearing loss.
- It is capable of providing more amplification as compared to other types.
- It is the best choice for hearing impaired children.
- It is larger in size hence more visible, but easier to handle, is sturdy and has a longer battery life.
- Requires batteries to be replaced from time to time.





B. RECHARGEABLE

- Suitable for all types of hearing loss.
- Long battery life for continuous use
- Fast charging.
- User-friendly design with easy controls and smartphone connectivity.
- Enhanced performance with noise reduction and wireless streaming.
- No battery replacements are required.















6. Rechargeable Earbuds

- Stylish Hearing Aids that don't look like Hearing Aids.
- Similar looking to Apple Air-pods or similar Bluetooth devices.
- Suitable for mild to severe hearing loss.
- Easy rechargeability with upto 26 hours use per charge.
- Portable charging case.
- Direct streaming from iPhone and iPad or android smartphones.



For Competitive Prices and Special Offers on **Hearing Aids**

Contact Centre for Hearing



www.centreforhearing.org



9811227269

Financing Options

Centre for Hearing has associations with financing companies offering 6-9-11-18 month 0% EMIs.